

“泰爾茂” 亞加尼針頭

“Terumo” AGANI Needle

衛署醫器陸輸字第000535號

使用前請務必詳閱原廠之使用說明書並遵照指示使用。

用途

本產品用途為與針筒或輸液套一起使用於藥品輸注或血液採集。

注意事項

無熱原產品，如包裝破損請勿使用。請勿儲存於極端溫度或濕度下。單次使用後請安全丟棄以避免感染風險。

產品規格

Code	Gauge	Length	Color	Bevel
AN*1838R	18G(1.2mm)	38mm	pink	Regular
AN*1938R	19G(1.1mm)	38mm	creamy	Regular
AN*2025R	20G(0.9mm)	25mm	yellow	Regular
AN*2038R	20G(0.9mm)	38mm	yellow	Regular
AN*2125R	21G(0.8mm)	25mm	deep green	Regular
AN*2138R	21G(0.8mm)	38mm	deep green	Regular
AN*2225R	22G(0.7mm)	25mm	black	Regular
AN*2232R	22G(0.7mm)	32mm	black	Regular
AN*2238R	22G(0.7mm)	38mm	black	Regular
AN*2338R	23G(0.6mm)	38mm	deep blue	Regular
AN*2325R	23G(0.6mm)	25mm	deep blue	Regular
AN*2332R	23G(0.6mm)	32mm	deep blue	Regular
AN*2425R	24G(0.55mm)	25mm	purple	Regular
AN*2438R	24G(0.55mm)	38mm	purple	Regular
AN*2516R	25G(0.5mm)	16mm	orange	Regular
AN*2525R	25G(0.5mm)	25mm	orange	Regular
AN*2613R	26G(0.45mm)	13mm	brown	Regular
AN*2713R	27G(0.4mm)	13mm	grey	Regular
AN*3013R	30G(0.3mm)	13mm	yellow	Regular

製造業者名稱：Zhejiang Kindly Medical Devices Co.,Ltd

製造業者地址：No. 758, 5th Binhai Road, Binhai Industrial Park,

Longwan District 325025 Wenzhou, Zhejiang Province, China

醫療器材商名稱：台灣泰爾茂醫療產品股份有限公司

醫療器材商地址：臺北市中正區忠孝西路一段4號8樓